

The gold standard of accuracy in the lab industry

The Vibrant Wellness **Wheat Zoomer™** uses proprietary technology to detect wheat sensitivity with >99% accuracy and 100% specificity, IgE, IgA, and IgG antibodies, celiac markers, and intestinal permeability.

Wheat Zoomer Markers List

- All known deamidated gliadins
- Alpha, alpha-beta, gamma and omega gliadin
- HMW and LMW glutenin family
- Zonulin protein
- Anti-zonulin
- Anti-actin
- Anti-lipopolysaccharide (LPS)
- tTG-DGP Fusion Peptides
- Wheat germ agglutinin (WGA)
- Differential transglutaminase: 2, 3 and 6
- Wheat IgE antibodies (for wheat allergies)
- Non-gluten wheat proteins
 - Farinins
 - Globulins
 - Serpins
 - Amylase/Protease Inhibitors

Published 11/18/2022